

ZIP kinase antibody

Catalog Number: 117059

Product name

ZIP kinase antibody

Specificity

Human; other species not tested.

Antibody description

ZIP kinase Rabbit Polyclonal antibody. Positive IF detected in HeLa cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001348).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using Catalog No:117059(DAPK3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)