

FAM111B antibody



Catalog Number: 110459

Product name

FAM111B antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human; other species not tested.

Validations

Antibody description

FAM111B Rabbit Polyclonal antibody. Positive IHC detected in human colon tissue, human prostate cancer tissue. Positive IF detected in HeLa cells. Positive WB detected in HeLa cells, human testis tissue, MCF7 cells. Observed molecular weight by Western-blot: 135kd




Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_198947). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110459(FAM111B antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human colon slide using Catalog No:110459(FAM111B Antibody) at dilution of 1:50

Ig Type

Rabbit IgG



Applications

ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

Immunohistochemistry of paraffin-embedded

FAM111B antibody

Catalog Number: 110459

human colon slide using Catalog
No:110459(FAM111B Antibody) at dilution of 1:50



Immunofluorescent analysis of Hela cells, using
FAM111B antibody Catalog No:110459 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit
IgG (red).