

HRASLS5 antibody



Catalog Number: 111459

Product name

HRASLS5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HRASLS5 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human heart tissue, human ovary tissue, human skin tissue, human spleen tissue. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of HRASLS5 recombinant protein (Accession Number: NM_001146729). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Immunohistochemical of paraffin-embedded human placenta using Catalog No:111459(HRASLS5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:111459(HRASLS5 antibody) at dilution of 1:50 (under 40x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111459(HRASLS5 Antibody) at dilution of

HRASLS5 antibody



Catalog Number: 111459

1:1000