

Anti-CAMKIV (3C10) Mouse antibody



Catalog Number: 168459

Product name

Anti-CAMKIV (3C10) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human

Applications

WB

Antibody description

CAMKIV (3C10) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human CaMKIV protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of CaMKIV in Romas and Jurkat cell lysates using CaMKIV mouse mAb (1:1000 diluted). Predicted band size: 52KDa. Observed band size: 55KDa.