

# SOD3 Polyclonal Antibody



Catalog Number: 161859

## Product name

SOD3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SOD3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-160 of human SOD3 (NP\_003093.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Immunofluorescence - SOD3 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using SOD3 antibody . Blue: DAPI for nuclear staining.



Immunofluorescence - SOD3 Polyclonal Antibody

Immunofluorescence analysis of PC-12 cells using SOD3 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.