

RTKN antibody



Catalog Number: 114859

Product name

RTKN antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RTKN Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in A431 cells, mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of RTKN recombinant protein (Accession Number: XM_017004635). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114859(RTKN antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:114859(RTKN antibody) at dilution of 1:100 (under 10x lens)