

Anti-ARHGEF10 antibody



Catalog Number: 175859

Product name

Anti-ARHGEF10 antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to ARHGEF10

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human arhgef10. aa range:339-388

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using ARHGEF10 Polyclonal Antibody