

Anti-CT A Rabbit antibody



Catalog Number: 166859

Product name

Anti-CT A Rabbit antibody

Specificity

Human

Antibody description

CT A Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CT A

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000-1/5000

Validations



Western blot detection of CT A in MCF7, HeLa cell lysates using CT A Rabbit pAb(1:1000 diluted). Predicted band size: 42KDa. Observed band size: 42KDa.