

MOAP1 Polyclonal Antibody



Catalog Number: 161159

Product name

MOAP1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to MOAP1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-351 of human MOAP1 (NP_071434.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:100

Validations



Immunohistochemistry - MOAP1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat kidney using MOAP1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - MOAP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using MOAP1 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - MOAP1 Polyclonal Antibody

MOAP1 Polyclonal Antibody



Catalog Number: 161159

Immunohistochemistry of paraffin-embedded

mouse brain using MOAP1 antibody at dilution of 1:100 (40x lens).