

ZNF311 antibody

Catalog Number: 117159

Product name

ZNF311 antibody

Specificity

Human; other species not tested.

Antibody description

ZNF311 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Positive IHC detected in human testis tissue, human kidney tissue. Observed molecular weight by Western-blot: 60-70 kDa

Preparation

This antibody was obtained by immunization of ZNF311 recombinant protein (Accession Number: NM_001382360). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:117159(ZNF311 antibody) at dilution of 1:4000



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:117159(ZNF311 Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human kidney slide using Catalog No:117159(ZNF311 Antibody) at dilution of 1:50