

# Anti-Cystatin C (6A4) Mouse antibody



Catalog Number: 168159

## Product name

Anti-Cystatin C (6A4) Mouse antibody

## Specificity

Human

## Antibody description

Cystatin C (6A4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant protein expressed in E.coli.

## Formulation

PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, IHC-P, ELISA

## Dilutions


WB: 1/1,000-2,000

IHC: 1/100-200 ELISA 1/10,000-20,000

## Validations



Western blot detection of Cystatin C in NCI-H460, SKOV-3, PC3 cell lysates using Cystatin C (6A4) Mouse mAb(1:1000 diluted). Predicted band size: 14KDa. Observed band size: 14KDa.



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.