

# ARNT2 antibody



Catalog Number: 108159

## Product name

ARNT2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ARNT2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Jurkat cells, Raji cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 80-85kd

## Preparation

This antibody was obtained by immunization of ARNT2 recombinant protein (Accession Number: BC036099). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108159(ARNT2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:108159(ARNT2 antibody) at dilution of 1:100 (under 10x lens)