

CRIP1 antibody



Catalog Number: 109559

Product name

CRIP1 antibody

Specificity

Human; other species not tested.

Antibody description

CRIP1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, PC-3 cells. Observed molecular weight by Western-blot: 8-9 kDa

Preparation

This antibody was obtained by immunization of CRIP1 recombinant protein (Accession Number: NM_001311). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109559(CRIP1 Antibody) at dilution of 1:600