

Anti-ELOA3P antibody



Catalog Number: 175959

Product name

Anti-ELOA3P antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ELOA3P

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human eloa3. aa range:381-430

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1:; Western Blot analysis of COLO205 cells using Elongin A3 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).