

ZNF211 antibody



Catalog Number: 116959

Product name

ZNF211 antibody

Immunogen

[Human ZNF211 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

ZNF211 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of ZNF211 recombinant protein (Accession Number: BC110856). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116959(ZNF211 Antibody) at dilution of 1:600