

Anti-PU.1/Spi1 Rabbit antibody



Catalog Number: 167959

Product name

Anti-PU.1/Spi1 Rabbit antibody

Specificity

Mouse

Antibody description

PU.1/Spi1 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of mouse PU.1/Spi1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/500

ICC/IF: 1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of PU.1/Spi1 in CHO-K1,C6 cell lysates using PU.1/Spi1 Rabbit mAb(1:1000 diluted).Predicted band size:31kDa.Observed band size:31kDa.