

# ALDH1A1 antibody



Catalog Number: 107959

## Product name

ALDH1A1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ALDH1A1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human hepat°Cirrhosis tissue. Positive WB detected in human liver tissue, human adrenal gland tissue, mouse liver tissue. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of ALDH1A1 recombinant protein (Accession Number: NM\_000689). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107959(ALDH1A1 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human liver using Catalog No:107959(ALDH1A1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:107959(ALDH1A1 antibody) at dilution of 1:50 (under 40x lens)