

# Anti-Caspase-6 Rabbit antibody



Catalog Number: 168959

## Product name

Anti-Caspase-6 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Caspase-6 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human Caspase-6

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of Caspase-6 in Jurkat, C6 cell lysates using Caspase-6 Rabbit mAb(1:1000 diluted). Predicted band size: 33kDa. Observed band size: 11kDa(cleavage), 33kDa.