

Human PSG9/PSG11 (His tag) recombinant protein



Catalog Number: 518959

General Information

Gene Name Synonym

Pregnancy-specific beta-1-glycoprotein 9, PS-beta-G-9, PSBG-9, Pregnancy-specific glycoprotein 9, PS34, Pregnancy-specific beta-1 glycoprotein B, PS-beta-B, Pregnancy-specific beta-1-glycoprotein 11, PS-beta-G-11, PSBG-11, Pregnancy-specific glycoprotein 11, Pregnancy-specific glycoprotein 7, PSG7, PSG11

Protein Construction

A DNA sequence encoding the Human PSG9/PSG11 Glu35-Ser426 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PSG9/PSG11 (His tag) recombinant protein