

HAS3 antibody



Catalog Number: 111259

Product name

HAS3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HAS3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, mouse liver tissue. Observed molecular weight by Western-blot: 63kd

Preparation

This antibody was obtained by immunization of HAS3 recombinant protein (Accession Number: NM_138612). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111259(HAS3 antibody) at dilution of 1:1000