

SRD5A1 Polyclonal Antibody



Catalog Number: 162259

Product name

SRD5A1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to SRD5A1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human SRD5A1 (NP_001038.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:100 - 1:200

Validations



Western blot - SRD5A1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SRD5A1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.