

# Kv4.2 antibody

Catalog Number: 112259

---

## Product name

Kv4.2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Kv4.2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, DU 145 cells, HeLa cells. Observed molecular weight by Western-blot: 70-80 kDa

## Preparation

This antibody was obtained by immunization of Kv4.2 recombinant protein (Accession Number: NM\_012281). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112259(KCND2 antibody) at dilution of 1:300