

# nNOS antibody

Catalog Number: 113259

---

## Product name

nNOS antibody

## Specificity

Human; other species not tested.

## Antibody description

nNOS Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000620).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113259(nNOS antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:113259(nNOS antibody) at dilution of 1:100 (under 40x lens)