

PTBP1 antibody



Catalog Number: 114259

Product name

PTBP1 antibody

WB: 1:1000-1:10000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

Antibody description

PTBP1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, BxPC-3 cells, HeLa cells, human lung tissue, K-562 cells, mouse lung tissue, rat heart tissue. Positive IHC detected in human prostate cancer tissue, human pancreas tissue. Positive IF detected in A2780 cells. Observed molecular weight by Western-blot: 59.6kd

Validations




Preparation

This antibody was obtained by immunization of PTBP1 recombinant protein (Accession Number: NM_002819). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:114259(PTBP1 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human prostate cancer slide using Catalog No:114259(PTBP1 Antibody) at dilution of 1:50

Ig Type


Rabbit IgG

Applications

ELISA, WB, IF, IHC

Dilutions

Recommended Dilution:



PTBP1 antibody



Catalog Number: 114259

Immunofluorescent analysis of A2780 cells, using PTBP1 antibody Catalog No:114259 at 1:25

dilution and Rhodamine-labeled goat anti-rabbit IgG (red).