

# SH3GL2 antibody



Catalog Number: 115259

## Product name

SH3GL2 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Immunogen

IP: 1:500-1:5000

[Human SH3GL2 Recombinant protein \(GST tag\)](#)

IHC: 1:20-1:200


## Specificity

Human, Mouse, Rat; other species not tested.

## Validations

## Antibody description

SH3GL2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 40 kDa




## Preparation

This antibody was obtained by immunization of SH3GL2 recombinant protein (Accession Number: BC032825). Purification method: Antigen affinity purified.

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115259(SH3GL2 antibody) at dilution of 1:50 (under 10x lens)

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115259(SH3GL2 Antibody) at dilution of

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IP, WB

## Dilutions

# SH3GL2 antibody



Catalog Number: 115259

IP Result of anti-SH3GL2 (IP:Catalog No:115259, 4ug; Detection:Catalog No:115259 1:1000) with mouse brain tissue lysate 4000ug.

