

SH3GL2 antibody



Catalog Number: 115259

Product name

SH3GL2 antibody

Immunogen

[Human SH3GL2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SH3GL2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of SH3GL2 recombinant protein (Accession Number: BC032825). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:115259(SH3GL2 antibody) at dilution of 1:50 (under 10x lens)



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115259(SH3GL2 Antibody) at dilution of

SH3GL2 antibody

Catalog Number: 115259

IP Result of anti-SH3GL2 (IP:Catalog No:115259, 4ug; Detection:Catalog No:115259 1:1000) with mouse brain tissue lysate 4000ug.

