

# TP53INP1 antibody

Catalog Number: 116259

## Product name

TP53INP1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TP53INP1 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human breast cancer tissue. Positive IF detected in HepG2 cells.

## Preparation

This antibody was obtained by immunization of TP53INP1 recombinant protein (Accession Number: NM\_033285). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF


## Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:20-1:200


## Validations



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:116259(TP53INP1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:116259(TP53INP1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells using Catalog No:116259(TP53INP1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)