

# Centrin 1 antibody



Catalog Number: 109259

## Product name

Centrin 1 antibody

## Immunogen

[Human CETN1 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

Centrin 1 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in HL-60 cells, human heart tissue. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of Centrin 1 recombinant protein (Accession Number: NM\_004066). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109259(CETN1 antibody) at dilution of 1:300



Immunofluorescent analysis of Hela cells, using CETN1 antibody Catalog No:109259 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).