

# Anti-nusA antibody



Catalog Number: 176659

## Product name

Anti-nusA antibody

## Specificity

E.coli

## Antibody description

Mouse monoclonal antibody to nusA

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at -20°C

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, ELISA

## Dilutions

WB: 1:2,000-1:5,000

## Validations




Fig1: Western blot analysis on E.coli cell lysates using anti- NusA Mouse mAb (Cat. # 176659#).