

Anti-nusA antibody



Catalog Number: 176659

Product name

Anti-nusA antibody

Specificity

E.coli

Antibody description

Mouse monoclonal antibody to nusA

Preparation

This antigen of this antibody was recombinant protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at -20°C

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, ELISA

Dilutions

WB: 1:2,000-1:5,000

Validations




Fig1: Western blot analysis on E.coli cell lysates using anti- NusA Mouse mAb (Cat. # 176659#).