

Anti-HOXB4 Rabbit antibody



Catalog Number: 168659

Product name

Anti-HOXB4 Rabbit antibody

Specificity

Human

Antibody description

HOXB4 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human HOXB4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/20-1/50

FC: 1/20

IP: 1/20

Validations



Western blot detection of HOXB4 in K562 cell lysates using HOXB4 Rabbit mAb(1:1000 diluted).Predicted band size:28kDa.Observed band size:34kDa.