

# C14orf156 antibody



Catalog Number: 108659

## Product name

C14orf156 antibody

## Specificity

Human; other species not tested.

## Antibody description

C14orf156 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 12 kDa

## Preparation

This antibody was obtained by immunization of C14orf156 recombinant protein (Accession Number: NM\_031210). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108659(C14orf156 Antibody) at dilution of 1:600