

HIST1H1T antibody



Catalog Number: 111359

Product name

HIST1H1T antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

HIST1H1T Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of HIST1H1T recombinant protein (Accession Number: NM_005323). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111359(HIST1H1T Antibody) at dilution of



1:300

mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:111359(HIST1H1T Antibody) at dilution of 1:300