

NUDT2 antibody



Catalog Number: 113359

Product name

NUDT2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NUDT2 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HEK-293 cells, mouse spleen tissue, mouse colon tissue. Positive IHC detected in human breast cancer. Observed molecular weight by Western-blot: 17kd

Preparation

This antibody was obtained by immunization of NUDT2 recombinant protein (Accession Number: NM_001161). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113359(NUDT2 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:113359(NUDT2 antibody) at dilution of 1:100 (under 25x lens)