

UGT1A1 Polyclonal Antibody



Catalog Number: 164759

Product name

UGT1A1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to UGT1A1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human UGT1A1 (NP_000454.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

Validations



Immunofluorescence - UGT1A1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using UGT1A1 antibody . Blue: DAPI for nuclear staining.