

RNF4 antibody



Catalog Number: 114759

Product name

RNF4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RNF4 Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of RNF4 recombinant protein (Accession Number: NM_001185009). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-RNF4 (IP:Catalog No:114759, 4ug; Detection:Catalog No:114759 1:300) with PC-3 cells lysate 3200ug.