

CACNA1G antibody



Catalog Number: 108759

Product name

CACNA1G antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CACNA1G Rabbit Polyclonal antibody. Positive IF detected in mouse heart tissue.

Preparation

This antibody was obtained by immunization of CACNA1G recombinant protein (Accession Number: NM_198382). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF

Dilutions


Recommended Dilution:

IF: 1:20-1:200

Validations

CACNA1G Polyclonal Antibody

Immunofluorescence of frozen mouse heart section. 1:100 dilution overnight at 4°C in Block. Secondary antibody: Alexa 594 Goat anti-rabbit. Apparent sarcolemmal expression, as expected.



IF result of anti-CACNA1G (Catalog No: 108759) in frozen mouse heart section by Dr Nicola Smart.