

DCTN4 antibody

Catalog Number: 109759

Product name

DCTN4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCTN4 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in human brain tissue, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 62 kDa

Preparation

This antibody was obtained by immunization of DCTN4 recombinant protein (Accession Number: NM_016221). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109759(DCTN4 antibody) at dilution of 1:800



IP Result of anti-DCTN4 (IP:Catalog No:109759, 4ug; Detection:Catalog No:109759 1:700) with Jurkat cells lysate 2800ug.