

TLE3 antibody

Catalog Number: 116069

Product name

TLE3 antibody

WB: 1:500-1:5000

IP: 1:200-1:1000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

Antibody description

TLE3 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A431 cells, HeLa cells, human lung tissue. Positive IP detected in HeLa cells. Positive IF detected in HepG2 cells. Positive IHC detected in human breast cancer tissue, human lung cancer tissue. Observed molecular weight by Western-blot: 83kd

Validations




Preparation

This antibody was obtained by immunization of TLE3 recombinant protein (Accession Number: BC015729). Purification method: Antigen affinity purified.

HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116069(TLE3 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, IP


Immunohistochemistry of paraffin-embedded human breast cancer slide using Catalog No:116069(TLE3 Antibody) at dilution of 1:50

Dilutions

Recommended Dilution:


TLE3 antibody

Catalog Number: 116069



Immunohistochemistry of paraffin-embedded human lung cancer slide using Catalog No:116069(TLE3 Antibody) at dilution of 1:50

IP Result of anti-TLE3 (IP:Catalog No:116069, 4ug; Detection:Catalog No:116069 1:300) with HeLa cells lysate 1600ug.



Immunofluorescent analysis of HepG2 cells using Catalog No:116069(TLE3 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG

