

Anti-APG5L/ATG5 Rabbit pA



Catalog Number: 167069

Product name

Anti-APG5L/ATG5 Rabbit pA

Specificity

Human, Mouse, Rat

Antibody description

APG5L/ATG5 Rabbit pA

Preparation

Antigen: A synthetic peptide of human APG5L/ATG5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000-1/2000

IHC: 1/20-1/50

ICC/IF: 1/20-1/100

FC: 1/50

Validations



Western blot detection of Atg5/Apg5L in Raji cell lysates using Atg5/Apg5L Rabbit pAb(1:1000 diluted).Predicted band size:33KDa.Observed band size:55KDa.