

Loncastuximab tesirine reference antibody



Catalog Number: 177069

Product name

Loncastuximab tesirine reference antibody

Specificity

Human CD19

Antibody description

Anti-CD19 Antibody (Loncastuximab tesirine)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications

ELISA

Validations



SDS-PAGE

Anti-CD19 Antibody (Loncastuximab tesirine) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.