

# HS2ST1 antibody

Catalog Number: 111469

## Product name

HS2ST1 antibody

## Specificity

Human; other species not tested.

## Antibody description

HS2ST1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in human liver tissue, human brain tissue, human heart tissue, PC-3 cells. Observed molecular weight by Western-blot: 27kd

## Preparation

This antibody was obtained by immunization of HS2ST1 recombinant protein (Accession Number: NM\_001134492). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111469(HS2ST1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:111469(HS2ST1 antibody) at dilution of 1:50 (under 10x lens)