

Phospho-Jak3-Y980/981 pAb



Catalog Number: 166469

Product name

Phospho-Jak3-Y980/981 pAb

Specificity

Human

Antibody description

Polyclonal antibody to Phospho-Jak3-Y980/981

Preparation

Antigen: A synthetic phosphorylated peptide around Y980 of human Jak3 (NP_000206.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-Jak3-Y980/981 pAb

Western blot analysis of extracts of 293T cells, using Phospho-Jak3-Y980/981 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.