

Anti-TP53INP2 antibody



Catalog Number: 176469

Product name

Rabbit IgG

Anti-TP53INP2 antibody

Applications

Specificity

WB, IHC-P

Human

Dilutions

Antibody description

WB: 1:500-1:2000

Rabbit polyclonal antibody to TP53INP2

IHC-P: 1:100-1:300

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human tp53inp2. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations




Fig1.; Western blot analysis of SH-SY5Y 293T mouse-kidney lysis using TP53INP2 antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).