

Anti-IKB alpha Rabbit antibody



Catalog Number: 167469

Product name

Anti-IKB alpha Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

IKB alpha Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human IKB alpha

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200


ICC/IF: 1/50-1/200

IP: 1/30

Validations




Western blot analysis of IKB alpha expression in



HeLa cell lysate.

Immunofluorescent analysis of HeLa cells, using IKB alpha Antibody .



Immunohistochemical analysis of paraffin-embedded mouse stomach, using IKB alpha Antibody.