

# phospho(403/404)-TDP43 antibody



Catalog Number: 107469

## Product name

Monoclonal

phospho(403/404)-TDP43 antibody

## Ig Type

Mouse IgG2a

## Specificity

Human; other species not tested.

## Applications

ELISA, WB

## Antibody description

phospho(403/404)-TDP43 Mouse Monoclonal antibody. Positive WB detected in K-562 cells, human placenta tissue. Observed molecular weight by Western-blot: 25 kDa

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007375).  
Purification method: Protein A purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107469(phospho(403/404)-TDP43 antibody) at dilution of 1:1000