

# phospho(403/404)-TDP43 antibody



Catalog Number: 107469

---

## Product name

phospho(403/404)-TDP43 antibody

## Specificity

Human; other species not tested.

## Antibody description

phospho(403/404)-TDP43 Mouse Monoclonal antibody. Positive WB detected in K-562 cells, human placenta tissue. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_007375). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107469(phospho(403/404)-TDP43 antibody) at dilution of 1:1000