

DLK1 Polyclonal Antibody



Catalog Number: 163869

Product name

DLK1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DLK1

Preparation

Antigen: Recombinant protein of human DLK1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:200 - 1:500

IHC 1:50 - 1:100

IF 1:50 - 1:100

Validations



Western blot - DLK1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using DLK1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.