

RUFY2 antibody



Catalog Number: 114869

Product name

RUFY2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RUFY2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, mouse lung tissue. Observed molecular weight by Western-blot: 65-70kd

Preparation

This antibody was obtained by immunization of RUFY2 recombinant protein (Accession Number: NM_001042417). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114869(RUFY2 antibody) at dilution of 1:1500