

# TBCA antibody

Catalog Number: 115869

## Product name

TBCA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TBCA Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells. Positive IHC detected in human prostate cancer tissue. Observed molecular weight by Western-blot: 13 kDa

## Preparation

This antibody was obtained by immunization of TBCA recombinant protein (Accession Number: NM\_004607). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115869(TBCA antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:115869(TBCA antibody) at dilution of 1:50 (under 10x lens)